

Genomic DNA from *Paenibacillus polymyxa*, Strain NCIB 8158

Catalog No. NR-52276

Product Description:

Genomic DNA was extracted from a preparation of *Paenibacillus polymyxa* (*P. polymyxa*), strain NCIB 8158. The bacterial preparation used for extraction of genomic DNA was produced by culture of BEI Resources NR-52263 lot 70033112. Genomic DNA was extracted using proprietary technology and is provided in TE buffer (10 mM Tris-HCl and 1 mM EDTA, pH 8).

Lot: 70033338

Manufacturing Date: 24MAR2020

| TEST | SPECIFICATIONS | RESULTS |
|--|--|---|
| Genotypic Analysis Digital DNA-DNA hybridization (dDDH) ¹ Sequencing of 16S ribosomal RNA gene (~ 1440 base pairs) | ≥ 70% for species identification ≥ 99% sequence identity to <i>P. polymyxa</i> , strain NCIB 8158 (GenBank: AFOX01000032.1) | <i>P. polymyxa</i> (99.6%) ^{2,3} <i>P. jamilae</i> (86.9%) ^{2,3} 99.4% sequence identity to <i>P. polymyxa</i> , strain NCIB 8158 (GenBank: AFOX01000032.1) |
| Agarose Gel Electrophoresis | High molecular weight chromosomal DNA | High molecular weight chromosomal DNA (Figure 1) |
| Concentration by PicoGreen® Measurement | 0.7 to 1.5 µg in 25 to 100 µL per vial | 1.2 µg in 33 µL per vial (37.5 µg per mL) |
| Amount per Vial | 0.7 to 1.5 µg | 1.2 µg |
| Functional Activity by PCR Amplification 16S ribosomal RNA gene | ~ 1500 base pair amplicon | ~ 1500 base pair amplicon |
| OD₂₆₀/OD₂₈₀ Ratio | 1.7 to 2.1 | 1.9 |
| Bacterial Inactivation 10% of total yield plated on agar for 14 days ^{4,5} | No viable bacteria detected | No viable bacteria detected |

¹Relatedness between bacterial strains has traditionally been determined using DDH. For additional information, refer to Auch, A. F., et al. "Digital DNA-DNA Hybridization for Microbial Species Delineation by Means of Genome-to-Genome Sequence Comparison." *Stand. Genomic Sci.* 2 (2010): 117-134. PubMed: 21304684. dDDH analysis was performed using the Type (Strain) Genome Server.


²The whole genome of *P. polymyxa*, strain NCIB 8158 was sequenced using the Illumina® MiSeq® system. *De novo* contig sequences were generated using Unicycler v0.4.8-beta.

³dDDH analysis of the type strains of *P. polymyxa* and *P. jamilae* shows an 86.9% sequence similarity, suggesting that *P. polymyxa* and *P. jamilae* are genomically the same species.

⁴14 days under propagation conditions

⁵An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative and Gram-positive bacteria.

Figure 1: Agarose Gel Electrophoresis



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder
Lane 2: ~ 200 ng of NR-52276

/Heather Couch/

Heather Couch

Program Manager or designee, ATCC Federal Solutions

21 JUL 2021

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

