

Mycobacterium bovis*, Strain 95-1315, Gamma-Irradiated Whole Cells*Catalog No. NR-59211**

This reagent is the tangible property of the U.S. Government.

Product Description:

Mycobacterium bovis, strain 95-1315 was grown to late-log phase in sodium pyruvate alanine salts medium and inactivated by exposure to a Cs source with the delivery of 1.44 to 2.4 mRads of ionizing radiation (dose range is provided due to variation across the sample chamber).

Lot: 70058662**Manufacturing Date: 30JAN2023**

Production and QC testing were performed by Colorado State University (CSU). The CSU documentation for pooled lot 23.M.bovis.95-1315.12.29.5.WCg and 23.M.bovis.95-1315.1.12.6.WCg is attached.

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



WORK SHEET FOR γ -IRRADIATED CELLS

General Information

Lot Number: 23.M.bovis.95-1315.12.29.5.WCg and 23.M.bovis.95-1315.1.12.6.WCg
Species: M. bovis
Strain: 95-1315

Growth Conditions: Medium: 7H9/OADC
of Fernbachs: 8 total
Inoculation Date: 12.29.22 Harvest Date: 1.12.23
Inoculation Date: 1.12.23 Harvest Date: 1.26.23
SOP # used for Harvest: PP003.3

Wet Weight: 48.4g

Date irradiated: 1.30.23 By: JB

Additional Notes: NA

Alamar Blue QC Information


Date Started: 2.14.23 Date Checked: 2.16.23

Positive Control: Positive X Negative


Sample: Positive Negative X

Aliquot Information

1 x 48.4g in 230ml falcon tube



(Research Associate) 2-16-23
(date)



(Laboratory Supervisor) 2/20/23
(date)

Initial 0 hr 570/600					
raw	2-14-23.	0.78	0.75	0.72	0.73
23.Mbovis.95-1315	48.4g	0.78	0.74	0.72	0.73
		0.77	0.73	0.72	0.72
		0.80	0.74	0.73	0.73

24hr 570/600					
	2-15-23.	0.81	0.73	0.73	0.73
		0.80	0.73	0.72	0.72
		0.79	0.73	0.73	0.73
		1.93	0.78	0.74	0.73

48 hr 570/600					
	2-16-23.	0.84	0.74	0.73	0.74
		0.84	0.75	0.73	0.73
		0.83	0.74	0.73	0.73
		3.78	0.83	0.75	0.74

Initial Minus Blank				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.06	0.03	0.01	0.01
	0.06	0.02	0.01	0.01
	0.05	0.01	0.01	0.01
Pos. Cont.	0.08	0.02	0.01	0.01

24hr Minus Blank				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.09	0.01	0.00	0.01
	0.07	0.01	0.00	0.00
	0.07	0.01	0.01	0.01
Pos. Cont.	1.21	0.06	0.01	0.01

48hr Minus Blank				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.11	0.01	0.00	0.01
	0.10	0.01	0.00	0.00
	0.10	0.01	0.00	0.00
Pos. Cont.	3.05	0.09	0.01	0.01

Initial Average				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.06	0.02	0.01	0.01
Pos. Ctrl.	0.08	0.02	0.01	0.01

24hr Average				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.08	0.01	0.00	0.00
Pos. Ctrl.	1.21	0.06	0.01	0.01

48hr Average				
	OD .1	OD .01	OD .001	OD.0001
Sample	0.10	0.01	0.00	0.00
Pos. Ctrl.	3.05	0.09	0.01	0.01

	Sample .1	Sample .01	Pos .1	Pos .01
0	0.06	0.02	0.08	0.02
24	0.08	0.01	1.21	0.06
48	0.10	0.01	3.05	0.09

